APPLICATION NOTE #1085

The Determination of Mercury in Soil and Water According to U.S. EPA SOW 846 &

Method 7473 by Thermal Decomposition, Amalgamation and Cold Vapor Atomic Absorption

INTRODUCTION

The accurate determination of mercury in soils and waters is critical to the remediation of contaminated sites. The *Hydra* II_c Direct Mercury Analyzer provides a simple and convenient way to analyze these materials in about 7 minutes without any sample digestion or creation of hazardous chemical waste to be disposed.

The **Hydra** II_c heats samples in an oxygen stream to decompose any mercury compounds which are then collected on an amalgamator. After all of the sample's mercury has been collected, the amalgamator is heated releasing the trapped mer-

cury for determination by cold vapor atomic absorption (Figure 1).



INSTRUMENTAL

Table 1 shows the instrument parameters employed for the calibration. Some of the sample matrices analyzed required slightly different drying and decomposition conditions for optimum results.





Table 1: System Parameters

Parameter	Water	Biological	Coal	Soil			
Dry	300°C for	300°C for	300°C for	300°C for			
Dry	30-140 sec.	30-140 sec. 70 sec.		30 sec.			
Decomposition	800°C for	800°C for	800°C for	800°C for			
	150 sec.	150 sec.	250 sec.	150 sec.			
Catalyst	600°C						
Catalyst Wait Peri-	60.000						
od	60 sec.						
Gold Trap	600°C for 30 sec.						
Measurement	100 sec.						
Oxygen Flow	350 ml/min						







SAMPLE PREPARATION AND CALIBRATION

The Hydra II_c employs two optical paths of differing lengths for extended dynamic range. Its software automatically selects the best analytical signal for each measurement. The calibration curves are displayed in Figures 3 and 4 as micro absorbance vs. total mercury injected.

These calibration curves were generated using weighed deposits of aqueous standards in concentrations of 0.1, 1.0 and 10.0 ppm (w/w) mercury. About 0.1–0.2 gm of each sample was placed into a tared nickel boat for analysis.

RESULTS


Before any analysis began the Instrument Detection Limit (IDL) and Method Detection Limit (MDL) were determined. The IDL was obtained by analyzing seven replicate blanks and the MDL was determined by analyzing seven replicates of NIST SRM8437 (Hard Red Spring Wheat Flour). Using 3.14 times the standard deviation*, the IDL was determined to be 0.005 ng and the MDL

was determined to be 0.02 ng. The data used to calculate the IDL and MDL are presented in Table 2. Certified reference materials (CRMs) for water, coal, biological tissue and

Sample Name Average sample Certified Concen-Measured Matrix weight tration Concentration (PPM) (PPM) (g) Water ERA WasteWatR 0.215 0.015 0.016 SABS SARM20 0.250 0.233 Coal 0.133 Biological NIST 1515 0.135 0.044 0.043 Soil NIST 2704 0.121 1.480 1.451

Table 3. Certified reference material results

soil matrices were analyzed. Certified and measured concentrations are provided in Table 3.



5.5765e-010

A =

Figure 4. Hydra II_c High Concentration Range (40–500 ng)

Replicate	IDL (ng Hg)	MDL (ng Hg)	
1	0.010	0.192	
2	0.011	0.179	
3	0.010	0.179	
4	0.014	0.189	
5	0.013	0.195	
6	0.011	0.180	
7	0.012	0.182	
Std Dev (S.D.)	0.0015	0.0067	
3.14(S.D.)	0.0047	0.021	

Table 2: IDL and MDL Results

*Student t-test for 99% distribution with seven replicates (6 degrees of freedom).

www.TeledyneLeemanLabs.com

Each sample type was analyzed in duplicate and a spike was analyzed to determine recovery. Sample spike (MS) recovery and spike duplicate (MSD) relative percent difference are pre-

sented in Table 4. Both the

spike recovery and the relative

percent difference between the spike and spike duplicate must be within ± 20%. The average sample weight analyzed also is displayed in Table 4.

Quality control checks were analyzed every ten samples for both the low and high curves. The method requires the results for the quality control checks to be within ±20%. Results were consistently within ±10%.

			А	В	С		
Sample Type	Sample Conc. (ng/g Hg)	Sample Weight (g)	Sample Mass (ng Hg)	MS Result (ng Hg)	Spike Added (ng Hg)	Recovery* (%)	Relative Percent Diff.** (%)
Waste	Waste 269 water	0.127	34.16	55.91	25	87.0	1.0
water		0.143	38.47	60.43	25	87.9	
Coal 27	27.1	0.138	4.0	33.63	31.75	94.1	7.4
	27.1	0.135	7.74	25.22	24.66	87.4	
Leaf	10.15	0.090	0.91	3.56	2.95	89.8	11.4
		0.081	0.82	5.41	4.57	100.4	
Soil	56.5	0.697	39.4	76.62	40.32	92.3	15.0
		0.727	41.1	88.5	44.1	107.5	

Table 4: Recovery and precision study

MS % Recovery is calculated by subtracting the sample result from the spike result and dividing the difference by the spike added. To report as a percentage, the calculated value is multiplied by 100. % Recovery = 100(B-A)/C.

**MSD Relative Percent Difference is calculated by multiplying the difference in the duplicate recoveries by 200 and dividing the result by the sum of the duplicate recoveries.

CONCLUSIONS

All the quality control requirements of U.S. EPA Method 7473¹ were satisfied. Sample results obtained showed excellent

precision and recovery. All CRM's showed good accuracy without any sample pretreatment and required only about 7 minutes per sample analysis. Proper drying and decomposition conditions produce a matrix independent analysis of mercury in many sample types.

MORE ABOUT THE HYDRA IIc

The **Hydra II**_c provides a 70 position autosampler and has onthe-fly loading capability for virtually unlimited capacity. Additionally, a conversion kit is available which employs chemical reduction to satisfy the monitoring of drinking water in accordance with USEPA Method 245.1 and European Standards EN1483 and EN13806.

¹www.epa.gov/osw/hazard/testmethods/sw846/pdfs/7473.pdf