



# RoHS

**EDXRF Bench-Top Spectrometer  
Pre-calibrated for RoHS**



  
**RoHS  
COMPLIANT**



Pre calibrated for restricted materials mandated under RoHS

Detection levels from ppm to 100%

Fast and non-destructive analytical method

Client server mode for remoted sites

# RoHS Spectrometer

The Restriction of Hazardous Substances (RoHS) is a set of standards that limit the use of certain toxic metals in electrical, electronic equipment and other consumables. System provides accurate analysis of restricted elements such as: Lead (Pb), Mercury (Hg), Cadmium (Cd), Chromium (Cr), Bromine (Br). The RoHS Spectrometer uses a high resolution detector, an integrated camera, and a powerful X-ray tube with variable spot sizes to accommodate samples of various sizes and to measure the existence of extremely low levels of restricted substances.

## Accurately analyze

Restricted elements Pb, Hg, Cd, Cr, Br

## Sample camera

Integrated camera and micro X-ray spot for full identification of the area of interest

## Ease-of-use

Automatic matrix identification allows operation of non-technical personnel with high degree of confidence

## RoHS designated software

Intuitive and user-friendly proprietary software

## Substances banned under RoHS

Lead (Pb), mercury (Hg), cadmium (Cd), hexavalent chromium (CrVI), polybrominated biphenyls (PBB), polybrominated diphenyl ethers (PBDE), and four different phthalates (DEHP, BBP, BBP, DIBP).

## RoHS standards reports

Producing accurate printed reports describing restricted elements detection level in reference to RoHS standards






# Easy-to-use Software

The designated software automatically detects the material type and selects the compatible testing methods. X-ray scan is being active and quickly a spectrum is presented and analysed. The tests are pre-calibrated according to RoHS standards and thresholds, while the resulting reports present the list of detected restricted elements concentration in compare to the specific defined range with accurate pass or fail results per element. A sample camera is constantly active, enabling to simultaneously view the examined sample spot location on the substrate. Each test can be individually saved into a library with a specific unique name, saving the spectrum and the analysis results.

This efficient pre-calibrated RoHS dedicated software provides an easy to use GUI and platform to operate and collect your RoHS tests.



## Key Applications

RoHS/WEEE compliance testing and screening of regulated elements (Pb, Hg, Cd, Cr, Br)





# RoHS Technical Specifications

	PD Version	SDD Version
<b>Detectable Range</b>	Pb, Hg, Cd, Cr, Br	Pb, Hg, Cd, Cr, Br [Na (11) - U (92)]
<b>Detectable Concentration</b>	ppm - 100%	
<b>X-ray Tube</b>	Mo - anode	
<b>X-ray Source</b>	50kV, 50W	
<b>Excitation Type</b>	Direct with filters	
<b>Spot Size</b>	Micro spot - Ø 1 mm. Macro spot - Ø 8 mm (On a sample)	
<b>Stability</b>	Precision 0.1% at ambient temperature	
<b>Detector</b>	PIN diode thermoelectrically cooled	SDD version
<b>Resolution (FWHM)</b>	155 eV ± 10eV at 5.9 keV.	135eV ± 5eV
<b>Autosampler</b>	One position	
<b>Work Environment</b>	Air/ Helium	
<b>Customized Tube Filters</b>	Six software selectable	
<b>Customized Collimator</b>	0.3 mm - 3 mm	
<b>Power Supply</b>	110-230VAC 50/60Hz	
<b>Pulse Processing</b>	Digital Pulse Processor	
<b>System Dimensions</b>	Unpacked: 55W x 55L x 32H, Packed: 80W x 80L x 65H	
<b>System Weight</b>	50kg (net), 90kg (gross)	
<b>Chamber Dimensions</b>	22 x 22cm, H=5cm	
<b>Camera</b>	CCD Camera	
<b>Computer</b>	Integrated PC	
<b>Operating Software</b>	User-friendly Operator mode (Microsoft Windows™ OS)	
<b>Control</b>	Automatic control of excitation, detection, sample handling and data processing	
<b>Spectrum Processing</b>	Automatic escape peak and background removal. Automatic peak deconvolution. Graphical statistics	
<b>Quantitative Analysis Algorithms</b>	Multi-element regression with inter-element corrections (six models available). Gross, net, fit and digital filter intensity methods	

## Distribuce v ČR a SR:

BAS Rudice s.r.o., U Vodárny 1, Blansko. ČR

Tel: +420 541 126 090

E-mail: bas@bas.cz

Xenometrix combines the latest technological developments with innovative engineering, to provide cost-effective solutions to a wide range of industries and applications.